



ORDERING INFORMATION

2082 - Intelect[®] RPW unit - 100 V 2073 - Intelect® RPW unit - 120 V 2074 - Intelect® RPW unit - 230 V

Standard Accessories:

28700 - Standard accessory kit that includes the following: 28699 - D-ACTOR® hand piece applicator 28736 - Projectile (x3) 17638 - R15 15 mm ESWT transmitter 28724 - D20-S D-ACTOR® 20 mm transmitter 28737 - Sealing set 28738 - Guide tube (x2) 28739 - Cleaning brush

4248 - Conductor™ transmission gel 250 ml (8.5 oz) bottle (x1) 27465 - Patient data card (x1) 28668 - User manual on CD

Optional Accessories:

28725 - D35 D-ACTOR® 35 mm transmitter 28726 - F15 Focus-Lens 15 mm transmitter 28728 - DI15 Deep Impact® 15 mm transmitter 28729 - C15 CERAma-x® ESWT 15 mm transmitter

> H F

28730 - V-ACTOR® hand piece applicator kit that includes the following: 28740 - V25 V-ACTOR® 25 mm vibration transmitter 28741 - V40 V-ACTOR® 40 mm vibration transmitter



www.medicaltronik.ca

DJO International 7 The Pines Business Park Broad Street Guildford Surrey GU3 3BH ENGLAND Tel: +44 (0)1483 459 659 Fax: +44 (0)1483 459 470

www.DJOglobal.eu





Shockwave by Chattanooga[™]





INTELECT® RPW

Intelect[®] RPW incorporates touch screen technology to ensure a high degree of simplicity. The easy-to-use menu-driven user interface guarantees reliable selection of all necessary parameters for treatment set-up as well as during patient treatment.

USER PROFILE

This Radial Pressure Wave clinical unit is ideally suited to the following applications:

- Orthopaedics
- Rheumatology
- Physical Medicine
- Physiotherapy
- Sport Clubs
- Acupuncture
- Aesthetic
- Veterinary

INDICATIONS

- Myofascial Trigger Points (MTrP)
- Localising and deactivating trigger points
- Activation of muscle and connective tissue
 - Increasing circulation
- Pulse vibration massage
- Disorder of tendon insertions
 - Plantar Fasciitis, heel pain, or heel spur
 - Tendinosis calcarea/supraspinatus-tendon
 - Radial and ulnar humeral epicondylitis
 - Achillodynia
 - Retropatellar pain syndrome
 - Tibial edge syndrome
 - Proximal iliotibial band friction syndrome/trochanteric
 - Insertional tendonitis

FEATURES & BENEFITS

1 2 CHANNEL OUTPUT CONNECTIONS

Intelect[®] RPW allows the connection of 2 hand pieces for quick and efficient use of complementary applicator transmitters during patient treatment. 2 hand piece applicators available:

D-ACTOR® Applicator

- **Applications:**
- Myofascial trigger points
- Other disorders of tendon insertions
- Activation of muscle and connective tissue
- Acupuncture shock wave therapy
- Cellulite treatment (stages I III)
- Muscle and connective tissue tightening
- Elimination of lymphatic congestion
- Body shaping

Medical massage

• Regulation of muscle tone



V-ACTOR® Applicator

The V-ACTOR® applicator is an ideal complement to Radial Wave Therapy and contributes to long-lasting treatment results.

Vibration treatment:

- Stimulates metabolism and accelerates the elimination of waste products
- Tightens tissue through muscular stimulation
- Improves skin firmness to restore a natural, younger appearance
- Offers a relaxing wellness experience

Applications:

- Body shaping
- Lymphatic drainage
- Wellness massage
- Muscle and connective tissue activation





2 STATE-OF-THE-ART NAVIGATION

This feature allows access to the desired functions at the touch of a finger. Treatment setup is quick and easy via Quick-link Indications, Clinical ProtocolsTM, User Protocols or manual adjustment:

• 200 user-defined protocol storage slots

• Over 80 clinical protocols on board



3 UNIQUE OPTIMAL ENERGY LEVEL ADJUSTMENT

- Unique ramp feature that allows gradual ramping of the intensity during treatment
- Provides an automatic and smooth intensity increase in tune with your patient
- Patient pain threshold can be recorded and treatment can be managed at the optimal painless energy level







4 CLINICAL RESOURCES

Intelect® RPW includes a full suite of reference resources available at the push of a button. These resources include:

- Full colour graphics library of anatomical images and pathologies
- Hand piece placement images • Radial pressure wave rationale
- Radial pressure wave contraindications







ADDITIONAL FEATURES & BENEFITS

• Convenient storage platform allows secure and easy access of hand pieces, transmitters and gel

 \leftarrow

- 360° swivel rotation of interface
- Large storage drawer for storage of clinical supplies
- Multilingual software
- Patient Documentation of treatment sessions, pre and post treatment patient pain scales, pain mapping and pain types all captured on patient data cards









