

# maxill

MYCOBACTERICIDAL • VIRUCIDAL† • BACTERICIDAL • FUNGICIDAL



# Minuteman NEX GEN DISINFECTANT WIPES



**tb Minuteman NEX GEN Disinfectant Wipes** are a new addition to maxill's line-up of hard surface disinfectants. When used as directed, tb Minuteman NEX GEN Disinfectant Wipes will disinfect in 1, 2 and 3 minutes against all the same microorganisms as tb Minuteman Hard Surface Disinfectant Towelettes. In addition to disinfecting against HIV-1, tb Minuteman Hard Surface Disinfectant Wipes have the added benefit of disinfecting against Hepatitis B (HBV) and Hepatitis C (HCV) to assure significant reduction of risk of cross-infection from these three commonly feared bloodborne viral pathogens in healthcare facilities and emergency vehicles.

The formulation for tb Minuteman NEX GEN Disinfectant Wipes has been designed using four active ingredients and supplemental ingredients to work synergistically to create an unparalleled broad-spectrum disinfectant wipe that is ideal for repeated daily use in dental offices, hospitals, nursing homes, and EMS & police vehicles. tb Minuteman NEX GEN Disinfectant Wipes may be used on a multitude of hard surfaces to inactivate a broad spectrum of pathogenic microorganisms.

- Broad Spectrum Bactericide, Virucide, Mycobactericide and Fungicide
- Meets CDC guidelines for intermediate level surface disinfection
- Cleaner and disinfectant in one product; save money, time and storage space
- Compatible with a wide range of hard, non-porous surfaces
- Low lint, heavy duty, non-woven, pre-saturated wipes, 6" x 7" and 7.5" x 9.75"
- Packaged in recyclable HDPE tubs with self-closing dispenser lids
- Convenient refill rolls are available to reduce waste and save you money

Designed with **You** in mind!

## Product Availability

160 wipes tub	#61160
160 wipes refill roll	#61190
110 jumbo wipes tub	#61180
110 jumbo wipes refill roll	#61195

Visit us online for more information, including SDS's and pricing:

[www.maxill.com/nexgen](http://www.maxill.com/nexgen)

### 1-minute contact time

Human Immunodeficiency virus type 1 (HIV-1)\*

### 2-minute contact time

Vancomycin Resistant *Enterococcus faecalis* (VRE)

Methicillin Resistant *Staphylococcus aureus* (MRSA)

Duck Hepatitis B Virus (Surrogate for Human Hepatitis B Virus)\*

Bovine Viral Diarrhea Virus (Surrogate for Human Hepatitis C Virus)\*

Avian Influenza A (H3N2) Reassortant virus

*Salmonella enterica*

### 3-minute contact time

*Mycobacterium bovis* BCG

*Pseudomonas aeruginosa*

*Staphylococcus aureus*

Poliovirus Type 1

*Trichophyton mentagrophytes*

Made in Canada with domestic and imported parts.

The information provided is given in good faith and believed to be accurate, but no warranty or guarantee is given or implied. Conditions of use under which the product is being applied will vary from the conditions of normal laboratory efficacy testing.

MP-REF-61160B  
Rev.0: 2019-11-15

maxill inc., St. Thomas, ON N5R 6C8  
[www.maxill.com](http://www.maxill.com) 1-800-268-8633